



For Tracheostomy Tubes


Speaking Valves



Heat/moisture exchanger (HMEs)




Occlusion Cap




NEW! AVAILABLE FROM 08/2014

Possible Combinations

Tubes, fenestrated



- **TRACOE twist and twist plus tubes, fenestrated, with 15 mm connector**
REF 302, REF 302-P, REF 304, REF 888-306, REF 888-306-P, REF 312, REF 312-P, REF 314, REF 888-316, REF 888-316-P




- **TRACOE comfort tube, fenestrated, with 15 mm connector**
REF 115, REF 215


Each in combination with:

- ◀ **phon assist I** (REF 650-T)
- ◀ **phon assist I with oxygen supply port** (REF 650-TO)
- ◀ **humid assist I** (REF 640-CT)
- ◀ **humid assist III** (REF 641)
- ◀ **humid assist IV** (REF 642)
- ◀ **phon assist II** (REF 655-T)
- ◀ **phon assist II (REF 655-T) + humid assist II** (REF 645)
- ◀ **Occlusion Cap** (REF 516)
- ◀ **22 mm adapter** (REF 622)


Tubes



- **TRACOE twist and twist plus tubes, with 15 mm connector**
REF 301, REF 301-P, REF 303, REF 305, REF 306, REF 309, REF 311, REF 311-P, REF 313, REF 316



- **TRACOE vario tubes, with 15 mm connector**
REF 450, REF 450-P, REF 451, REF 451-P, REF 455, REF 460, REF 460-P, REF 461, REF 461-P, REF 470, REF 470-P, REF 471, REF 471-P



- **TRACOE comfort tubes, with 15 mm connector**
REF 103-A, REF 104-A, REF 114, REF 105, REF 105-D, REF 204-A, REF 205, REF 214


Each in combination with:

- ◀ **humid assist I** (REF 640-CT)
- ◀ **humid assist III** (REF 641)
- ◀ **humid assist IV** (REF 642)
- ◀ **22 mm adapter** (REF 622)



Caution: All tubes with cuff have to be deflated before using a speaking valve!

Tubes, for neonates and children



- **TRACOE mini tubes with 15 mm connector**



REF 350,




REF 355

Each in combination with:


- ◀ **humid assist mini** (REF 643)
- ◀ **22 mm adapter** (REF 622)

For TRACOE buttons, Silicone Short Tubes and Adhesive Carriers

Speaking Valves



Heat/moisture exchanger (HMEs)




Stopper




Possible Combinations

Carrier with 22 mm housing



Each in combination with:

- ◀ **phon assist I** (REF 650-S)
- ◀ **phon assist I with oxygen supply port** (REF 650-SO)
- ◀ **humid assist I** (REF 640-S)
- ◀ **phon assist II** (REF 655-S)
- ◀ **phon assist II** (REF 655-S) + **humid assist II** (REF 645)
- ◀ **Stopper** (REF 621)



Each in combination with:

- ◀ **humid assist I** (REF 640-S)

